

**UNMOUNTED GENESIS MED & LARGE STENTS (GREY BOX)**

Reference #	STENT LENGTH UNEXPANDED MM	EXPANDED STENT (DIAMETER) MM	STENT LENGTH EXPANDED	Introducer Size [Fr]	AVERAGE % STENT SHORTENING
<b>PG184B (med)</b>	<b>18</b>	<b>5</b>	<b>18</b>	<b>6</b>	<b>6%</b>
		<b>6</b>	<b>17</b>		
		<b>7</b>	<b>17</b>		
		<b>8</b>	<b>15</b>		
<b>PG244B (med)</b>	<b>24</b>	<b>5</b>	<b>24</b>	<b>6</b>	<b>4%</b>
		<b>6</b>	<b>23</b>		
		<b>7</b>	<b>23</b>		
		<b>8</b>	<b>21</b>		
<b>PG295B (large)</b>	<b>29</b>	<b>5</b>	<b>28</b>	<b>7</b>	<b>9%</b>
		<b>6</b>	<b>28</b>		
		<b>7</b>	<b>27</b>		
		<b>8</b>	<b>26</b>		
		<b>9</b>	<b>25</b>		
		<b>10</b>	<b>24</b>		
<b>PG595B (large)</b>	<b>59</b>	<b>5</b>	<b>59</b>	<b>7</b>	<b>5%</b>
		<b>6</b>	<b>58</b>		
		<b>7</b>	<b>57</b>		
		<b>8</b>	<b>56</b>		
		<b>9</b>	<b>53</b>		
		<b>10</b>	<b>59</b>		