

VALVE TABLE SET UP

Valve table Setup: Edwards

- 2 small/wide tables taped at legs, table cover
- Metz and pick up from valve wash bowl set
- 4 large bowls and 4 small bowls
- 2 stop cocks
- 60cc syringe x 1
- 30cc syringe x 2
- 3cc LL syringe x 1
- Labels and pen
- 1 pack of towels
- Gowns and gloves (valve nurse will assist in gowning and gloving)
- Turp tubing for saline transferred from scrub table
- 500cc of Hep-saline in one large valve bowl
- NaCl injectable 1L in large bowls
- Visopaque x 1 small bottle in small bowl (**)

Valve table Setup: Evolute R

- 2 small/wide tables taped at legs, table cover
- Metz and pick up from valve wash bowl set
- 1 20 cc syringe
- Labels and pen
- 1 pack of towels
- Gowns and gloves (valve nurse will assist in gowning and gloving)
- Turp tubing for saline transferred from scrub table
- 3 One liter bags of normal saline (move from fridge to freezer the morning of)
- One bag transplant NS prn
- Visopaque x 1 small bottle in small bowl (**)